





See 3 mm next page.

Maximum PROBE Speeds (kHz)	
DI 4	18
XC4	22
XC5	18
5 SS	18
5 HS	14
5 Std	9



A - Si₃N₄ Thick Wall | B - Si₃N₄ Thin Wall | C - Zr Thick Wall | D - Zr Thin Wall

Maximum PROBE Speeds (kHz)	
XC7	12
7 SS	12
7 HS	9
7 Std	6
XC10	8.5
10 SS	8.5



MAS Turbine Cap Spinning Speeds

Maximum Spinning Speeds (kHz) For Caps at Room Temperature

Cap Style	4 mm	5-mm XC or SuperSonic	5-mm Standard & High-Speed	7-mm XC or SuperSonic	7-mm Standard & High-Speed	10-mm XC, SuperSonic
Kel-F	11	10	9	7	6	5
Torlon or GFT	22	18	14	12	9	8.5
Vespel	21	16	14	12	9	8
Aurum	22	18	14	12	9	8.5
Caps with o-rings	-----	10	9	7	6	5
Vespel w/screw	-----	-----	9	12	11	8

Note: This chart represents only material characteristics for caps. Check the Probe Specifications. The spinning speed is often more limited by the probe or the rotor material.



3 mm Maximum MAS Spinning Speeds

Use the lower of the speeds listed: considering the rotor, the cap, the temperature and the density maximum speeds

Rotors: <i>The maximum speed must be reduced as the density of the sample increases</i>			
DI3 Silicon Nitride Rotors Maximum Speed	28 KHz	For sample density = 1	
	26 KHz	For sample density = 3	
Maximum Spinning Speeds (kHz) For Caps at Room Temperature			
Cap Style	DI3	Spinning Temperature	Cap Material Temperature Range
Torlon Front Turbine Cap	28 KHz	At Room Temperature	-30° to 80 °C
Torlon Rear Tip Cap	28 KHz		
GFT* Front Turbine Cap	26 KHz	At Room Temperature	-120° to 160 °C
GFT* Rear Tip Cap	26 KHz		
Aurum Front Turbine Cap	18 KHz	At Room Temperature	-30° to 80 °C
Aurum Rear Tip Cap	18 KHz		
Kel-F Front Turbine Cap	11 KHz	At Room Temperature	-20° to 70 °C
Kel-F Rear Tip Cap	11 KHz		

Maximum Spinning Speeds (KHz) For Caps at Extended Temperatures

Glass Filled Torlon (GFT)*	14 KHz	At -80°C	-120° to 160 °C
Glass Filled Torlon (GFT)*	14 KHz	At +120°C	-120° to 160 °C

* Note: **GFT caps** can be used up to 250 °C or down to -170 °C if they are glued in with epoxy. However the probe must be rated for these extended temperatures.