Low-E Inverse Detection MAS Probe and Gradient Spectroscopy in Solid State NMR.

Bibhuti B. Das¹, Leonard Mueller², Laura Holte¹, Daniel Arcos¹, Marc Bremmer¹, John Staab¹, J. B. Spitzmesser¹, Paul D. Ellis¹, and F. David Doty¹

Doty Scientific, 700 Clemson Rd. Columbia, SC 29229.² Department of Chemistry, University California, Riverside, CA.

Introduction

Α

¹H detected solid-state NMR experiments have been rapidly evolving in past decade enabling sensitivity enhanced experiments to be applied in solid samples. On the other hand pulse field gradients (PFGs) are routinely used in solution NMR facilitating rapid acquisition of multi-dimensional experiments with efficient coherence selection pathways, diffusion studies and in magnetic resonance imaging as well. By combining the above two techniques it is possible to broaden solid state NMR applications beyond its current capabilities. Hence, there has been an increasing focus at Doty Scientific on design and development of high-resolution probes for moderate and fast magic angle spinning with the inclusion of gradient coils. Here the design, construction and performance of at least two different solid state probes equipped with gradient coils: (1) Magic angle spinning-magic angle gradient (MAS-MAG) probe optimized for inverse 1H **detection**; (2) Static triple resonance probe for solid state NMR of mechanically and magnetically oriented samples. Further, these probes have sufficient tuning flexibility to work efficiently with high loss biological samples and less sample heating. Additionally, gradient coils producing fields up to 500 G/cm and fast recovery time enables acquisition of high-resolution ¹H spectra for combined rotation and multiple pulses (CRAMPS) experiments under moderate MAS. Experimental results obtained from a stationary probe operating at 900 MHz ¹H frequency and new methodologies for gradient enhanced spectroscopy are demonstrated.

DOTY Low-E MAG-MAS Probe



Coherence selection using gradients in MAS-NMR

Solution NMR experiments used for coherence selected pathways using PFGs are modified to accommodate solid samples under MAS. Amplitude and phase modulated Lee-Goldburg RF pulses are applied for homonuclear dipolar decoupling during windowed acquisition to compensate for phase transients and to improve spectral resolution.



- Spectral resolution of 0.015 ppm.

¹H detected magic angle setting

High accuracy magic angle setting with proton detection is obtained using a sample with carbon labeled chloroform dissolved in partially oriented liquid crystalline (LC) matrix. The director(n) of the LC matrix align along the magnetic field when stationary, but aligns along the axis of rotation under rotor spinning. Here, scaled dipolar couplings are measured to monitor rotor angles.

Stimulated Echo PFG- diffusion NMR under MAS

Liquid crystals are viscous semi solids. They are often used in commercial products, chemical synthesis and as orienting media for NMR study of small molecules. A combination of MAS and PFG techniques are applied to study self diffusion of liquid crystals and the solutes in nematic and in isotropic phases. The technique enables measurement of varying diffusion coefficients in solids with respect to molecular sizes, viscousity, and temperature, suggesting a wider application in polymer chemistry and biological sciences.



300 MHz ¹H signal of chloroform in nematic 5CB. A) ¹³C decoupled 10 ms spin echo signal obtained under MAS. B) 1ms echo signal under MAS and no ¹³C decoupling applied. C) ¹³C coupled 1 ms echo signal under stationary conditions. 3 kHz splitting in C accounts for sum of J_{CH} -coupling and twice d_{CH} -coupling.

Water Suppression using RF and PFGs

Efficient solvent suppression is achieved by applying long low power RF and gradient pulses. Spectra are shown for glycine dissolved in deuterated water. Residual water signal is noted at 5 ppm and glycine peak near 3 ppm.



A) Stimulated echo pulse sequence for PFG-diffusion NMR study of a nematic liquid crystal 5CB under MAS. B) Array of ¹H spectra recorded under 7 kHz MAS with varying gradient strength. Black and red colors indicate isotropic and nematic phase of 5CB. C) Graphical demonstration of signal decay for solvent (5CB) and the solute (CHCl₃).

Acknowledgement

NIH

DSI acknowledges NIGMS-NIH for research funding grants R44GM108147-03, R44GM079888, and their dedicated program officials. DSI also acknowledges NSF sponsored National High Magnetic Field Laboratory for allowing use of 900 MHz widebore magnet and research scientists, Dr. Zheong Gan, Dr. Ivan Hung, Dr. Peter Gorkov and Dr. Jason Kitchen, for their technical help.

A) Radio frequency pulse sequence used for ¹H detection and homo-nuclear chemical shift correlation. B) Gradient recovery time of the probe for a 2 ms gradient pulse (100 G/cm) applied. Note that the probe recovers within 50 micro sec. C and D) Gradient assisted double and single quantum filtered spectra of Glycine powder recorded under 7.5 kHz MAS on a 300 MHz spectrometer. E) Two-dimensional ¹H- ¹H chemical shift correlation spectrum of Glycine. Direct detected 1D spectra are shown on top and left.

Gradient selected echo/anti-echo phase cycling in triple-resonance experiments at 900 MHz ¹H NMR.





References

•W. E. Maas, A. Bielecki, M. Ziliox, F. H. Laukien, and D. G. Cory, Magnetic Field Gradients in Solid •State Magic Angle Spinning NMR, Journal of Magnetic Resonance 141, 29–33 (1999). •C.A. Fyfe, J. Skibsted, H. Grondey, H. Meyer zu Altenschildesche, Pulsed field gradient multiplequantum MAS NMR spectroscopy of half-integer spin quadrupolar nuclei, CPL, 281, 44-48 (1998).

Triple resonance pulse sequence (A) used for ¹³C edited ¹⁵N detection (B) and ¹³C-¹⁵N correlation (E) under gradient selected echo/anti-echo phase cycling. C) Echo and D) anti-echo (red). Sample: ¹³C, ¹⁵N labeled Valine powder.