

## Dear Colleague,

Many of you need NMR probes and MRI coils for Agilent systems. Doty Scientific has been making high performance NMR probes and MRI coils over 30 years that are compatible with Agilent, Bruker, JEOL, Tecmag, and other systems.

At conferences this spring and summer, we'll be discussing ways to make DNP affordable. We'll be sharing information on our CryoMAS-DNP probe; our new Static DNP probe (see below); and a wide range of NMR probes and MRI imaging coils. Let's talk about your interests.

David and Judy Doty

## Upcoming Conferences

SUN  
19  
APR  
2015

FRI  
24  
APR  
2015

### ENC


Asilomar Conference Center  
Pacific Grove, California, USA  
[www.enc-conference.org](http://www.enc-conference.org)  
Visit our ENC Hospitality Suite  
'Toyon'

SUN  
31  
MAY  
2015

SUN  
5  
JUNE  
2015

### ISMRRM

Metro Ontario Convention Center  
Toronto, Ontario, Canada  
[www.ismrm.org/15](http://www.ismrm.org/15)  
Visit our Booth # 708




CryoMAS™ DNP Probe

## MRI Coils

Small Animal Imaging Platforms

RF Imaging Modules

Surface Coils



## NMR Probes

CryoMAS DNP Probes

Static DNP Probes

Solids / Liquids MAS Probes

Magic Angle Gradient Probes

PFG / Diffusion Probes

Liquids HR Probes

High Temperature Probes

## Static-DNP Probe

Making it possible for NMR groups to bring static H/X/Y/e<sup>-</sup> DNP into their labs


For both solids and liquids at high fields

The estimated price including the solid-state  $\mu\text{w}$  source, the Doty DNP probe, waveguides, and transitions, for 7 T to 11.7 T, is anticipated to be in the range of \$180K to \$220K

Adequate solid-state  $\mu\text{w}$  sources are available up to 11.7 T

There has been a lot of excitement about MAS-DNP over the past four years (e.g. its potential for two orders of magnitude increase in S/N), but the cost has appeared to be well beyond what most NMR labs could begin to consider – due to the need for a special magnet with superconducting sweep coils and an expensive gyrotron. However, you will probably be surprised to learn that usually less than 0.1% of the microwave ( $\mu\text{w}$ ) power from the gyrotron ends up being dissipated in the sample by mechanisms that represent fundamental limitations. This suggests that in principle it should be possible to do DNP (at the same temperature, 90-150 K) with  $\sim 0.1\%$  of the amount of  $\mu\text{w}$  power that has been thought necessary. If true, gyrotrons would not be needed. Low-cost solid-state sources could be used, and their bandwidth is more than sufficient to eliminate the need for superconducting sweep coils and a new NMR magnet.

Acknowledgement: Research reported here was supported by NIA of the National Institutes of Health under award numbers 1R41AG050394. The content is solely the responsibility of the authors and does not necessarily represent the official views of the NIH.



H/X/Y/e<sup>-</sup> Static DNP Probe

We have been hard at work for two years toward this vision and have made amazing progress. Our latest (fully detailed) simulations are saying our first-generation cavity for static (non spinning) H/X/Y/e<sup>-</sup> DNP at 7T will be about two orders of magnitude more efficient than what has previously been seen in high-field DNP probes. This revolutionary advance (patent pending) will make it possible for all NMR labs to begin applying high field H/X/Y/e<sup>-</sup> static DNP to their solids-NMR problems, such as in structure-function studies of membrane proteins.

And what about the potential for similar efficiency breakthroughs in cavities compatible with MAS-DNP? They're coming too – within about two years. Come by our 'Toyon' ENC Hospitality Suite to learn more, and get in the queue for more than ice cream – the next big wave in NMR !