

Two-coil Ultra High Field No-E H/X/Y for XVT MAS and Wideline

Laura L Holte, F. David Doty, Dylan Harrell, JB Spitzmesser, Daniel Arcos, and Paul D Ellis
Doty Scientific, Inc., Columbia, SC USA



INTRO Obtaining useful NMR data on low-gamma quadrupolar nuclides in solid samples remains a major challenge that in many cases will only be addressed with the combination of higher B_0 fields, larger samples, lower temperatures, and efficient triple-resonance methods in standard-bore magnets. That need served as the primary motivation for this development effort. While single-sample-coil (single-coil) circuits have usually been used at fields above 600 MHz, such circuits are incompatible with high rf decoupling in large samples at high fields. Hence, a two-coil approach was chosen with an inner solenoid and an outer segmented transverse proton coil. Because of its half-turn capacitive segmenting, the High Frequency (HF) E field in the sample is much lower than that seen in most “E-Free” probes, and it vanishes in the center – hence, the name “No-E”.

METHODS

- The inner X/Y solenoid is secured to a Si_3N_4 ceramic stator for minimal tuning change with temperature and with long-term aging.
- A novel method for tuning the outer resonator from ^1H to ^{19}F , suitable for all fields and sample sizes, was implemented and modeled.
- Accommodating the full range of nuclides, sample sizes, and polarizing fields with high efficiency and power handling requires lower stray capacitance and more flexibility in plug-in component options than appears possible with plug-in-circuit-board approaches.
- Individual plug-in tuning elements work better!
- Detailed rf circuit models for the solenoid X/Y circuit and the HF circuit were developed, based on an earlier similar structure, with tuning data from several sample sizes at fields from 7-14 T.
- The detailed models showed that with several changes in the probe structure and circuits, efficient fully multi-nuclear H/X/Y performance could be expected to at least 30 T with small or large samples.




Fig. 1. Highly simplified low-field approximation of the basic X/Y solenoid rf circuit for conceptual purposes, with approximate values for the 5-mm 9.4-T $^{31}\text{P}/^{13}\text{C}$ case.

- For best performance, the high-current X/Y plug-ins must be just below the sample region. They provide voltage/impedance reduction to the long leads going to the LF/MF variables farther down.
- The HF fine-tune variables are in the region below the main LF/MF plug-ins. Long leads to the HF resonator aren't a problem, as lead VSWR is kept sufficiently low for low lead losses.
- The LF/MF tune/match variable (quartz) capacitors are in the lowest rf zone.
- With ten plug-in capacitors and coils, detailed circuit simulations show the novel X/Y circuit achieves optimal performance for H/X and H/X/Y combinations from ^{31}P to below ^{107}Ag .

