



## Wideline Probe Accessories (1998-Present)

(Please provide probe serial number to be sure the correct parts are supplied)

(Prices listed are USA prices. Foreign prices higher, plus taxes.)

5mm (When ordering, please provide probe serial number.)										
Wideline Sample Containers			* WL End Caps For Glass Containers				WL End Caps For Plastic & Ceramic Containers			
Part #	Material	Price	Part #	Material	Type	Price	Part #	Material	Type	Price
5152	Glass *	\$ 45	5832	Kel-F	Short	\$ 30	5834	Kel-F	Short	\$ 30
5153	Zirconia	160	5833	Kel-F	Long	30	5835	Kel-F	Long	30
5154	Kel-F	60	5836	Aurum	Short	40	5838	Aurum	Short	40
5155	Aurum	60	5837	Aurum	Long	40	5839	Aurum	Long	40
5156	Peek	60	5841	Peek	Short	30	5843	Peek	Short	30
-	-	-	5842	Peek	Long	30	5844	Peek	Long	30

7 mm (When ordering, please provide probe serial number.)										
Wideline Sample Containers			* WL End Caps For Glass Containers				WL End Caps For Plastic & Ceramic Containers			
Part #	Material	Price	Part #	Material	Type	Price	Part #	Material	Type	Price
8256	Glass *	\$ 45	8116	Kel-F	Short	\$ 30	8118	Kel-F	Short	\$ 30
8257	Zirconia	160	8117	Kel-F	Long	30	8119	Kel-F	Long	30
8258	Kel-F	60	8132	Aurum	Short	40	8134	Aurum	Short	40
8259	Aurum	60	8133	Aurum	Long	40	8135	Aurum	Long	40
8260	Peek	60	8136	Peek	Short	30	8138	Peek	Short	30
-	-	-	8137	Peek	Long	30	8139	Peek	Long	30

8 mm (When ordering, please provide probe serial number.)										
Wideline Sample Containers			* WL End Caps For Glass Containers				WL End Caps For Plastic & Ceramic Containers			
Part #	Material	Price	Part #	Material	Type	Price	Part #	Material	Type	Price
8251	Glass *	\$ 45	8103	Kel-F	Short	\$ 30	8105	Kel-F	Short	\$ 30
8252	Zirconia	160	8104	Kel-F	Long	30	8106	Kel-F	Long	30
8253	Kel-F	60	8107	Aurum	Short	40	8109	Aurum	Short	40
8254	Aurum	60	8108	Aurum	Long	40	8111	Aurum	Long	40
8255	Peek	60	8112	Peek	Short	30	8114	Peek	Short	30
-	-	-	8113	Peek	Long	30	8115	Peek	Long	30

10 mm (When ordering, please provide probe serial number.)										
Wideline Sample Containers			* WL End Caps For Glass Containers				WL End Caps For Plastic & Ceramic Containers			
Part #	Material	Price	Part #	Material	Type	Price	Part #	Material	Type	Price
9501	Glass *	\$ 45	9251	Kel-F	Short	\$ 30	9253	Kel-F	Short	\$ 30
9502	Zirconia	160	9252	Kel-F	Long	30	9254	Kel-F	Long	30
9503	Kel-F	60	9255	Aurum	Short	40	9257	Aurum	Short	40
9504	Aurum	60	9256	Aurum	Long	40	9258	Aurum	Long	40
9505	Peek	60	9259	Peek	Short	30	9262	Peek	Short	30
-	-	-	9261	Peek	Long	30	9263	Peek	Long	30

\* The end caps for glass are different than for plastic and ceramic. Be sure to use the correct caps, otherwise it is likely that the glass containers will be broken.

Previous Generation Wideline Probe Accessories (prior to 1998) (When ordering, please provide probe serial number.) (Prices listed are USA prices. Foreign prices higher, plus taxes.)					
5 mm Part #	7 mm Part #	8 mm Part #	10 mm Part #	Item	Price Each
5151	7512	8250	9500	Pyrex or Ceramic Sample Container (5-mm and 7-mm are ceramic)	\$ 45
5720	7720	8110	9250	Short Kel-F plug cap	35
5760	7760	8100	9260	Long Kel-F plug cap	35
-	-	8170	9290	Short Teflon plug cap	35
-	-	8160	9300	Long Teflon plug cap	35
5870	-	8101	9400	Macor plug cap	60
5780	7780	8102	9350	Vespel or Aurum plug cap - short	40
5820	7820	-	-	Vespel or Aurum plug cap - long	40
96220	96240	96240	96174	Puller for Plug Caps	160

**If unsure about correct supplies, contact us with the DSI-serial number and we can help.**

Sales: [sales@dotynmr.com](mailto:sales@dotynmr.com), Sales: 803-738-6832, Main Office: 803-788-6497