









General Accessories (Air Connectors, Paints, Glues)

Item			Part #	Price US\$
 <p># 96383 Brass Snap Tite Probe Air Connector with male quick connect</p>	<p>#96383-XVT This XVT version is used for low temperature experiments – without a Doty cold gas VT system. Please supply the DotySerial number “DSI-__” when ordering.</p>			\$ 140
			# 96383-XVT	(XVT) \$ 165
<p style="text-align: center;">Filtered Air Line with Female Quick Connect These airlines are used for room temperature air cooling and other non-MAS probe air requirements.</p>  <p style="text-align: center;"># 96382</p>			# 96382	\$ 82
<p style="text-align: center;">Filtered Air Line with Female Quick Connect – For MAS These MAS air lines have internal channels that reduce turbulence at the drive and bearing inlets on MAS probes.</p>  <p style="text-align: center;"># 96382-MAS</p>			# 96382-MAS	\$ 121
<p style="text-align: center;">Filtered Air Line with Male Quick Connect</p>  <p style="text-align: center;"># 90617</p>			# 90617	\$ 82
<p style="text-align: center;">Quick Connect Set</p>  <p style="text-align: center;"># 96390</p>	<p style="text-align: center;">Female Quick Connect</p>  <p style="text-align: center;"># 58560</p>	<p style="text-align: center;">Male Quick Connect</p>  <p style="text-align: center;"># 58550</p>	# 96390	\$ 66
			# 58560	\$ 33
			# 58550	\$ 33
<p style="text-align: center;">Ball and Clamp Probe Air Connector to male connect (Used on Early Probes)</p> 			# 96381	\$ 82
Optical <u>Rotor</u> Paints			# 96109	\$ 95
Glue for Caps for XVT			# 99529	\$ 66

(US\$ – Foreign prices higher, plus taxes.)